IT

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
Levenhuk, Inc. 924-D East 124th Ave. Tampa, FL 33612 USA
Levenhuk® is a registered trademark of Levenhuk, Inc.
© 2006–2018 Levenhuk, Inc. All rights reserved.
www.levenhuk.com

20180618

levenhuk
Zoom&Joy

EN

The Levenhuk MS3 mechanical scale is designed to move microscope slides on the stage along X axis (left-right) and Y axis (forwards-backwards). It also helps determine slide coordinates. Each coordinate scale is accompanied by an auxiliary nonius scale, which allows for high precision measurements. The scale is recommended for use with metal stages; plastic stages are not sturdy enough to hold the scale and could be damaged.

Installation: screw the slide holders out of the standard stage using a screwdriver. Install the measuring scale on the stage and fix it with the stage screw (5).

The scale is compatible with Levenhuk S NG, Rainbow PLUS microscopes. The installation is not recommended for microscopes that have plastic bodies as damage to the microscope may occur.

Movement range along X axis, mm	0–57
Movement range along Y axis, mm	0–30
Measuring scale division value, mm	1
Distance between the dowel pins, mm	35±0.1
Distance between the dowel pins and the stage screw, mm	17.5±0.1
Dowel pin diameter, mm	3
Stage screw thread size	M4
Scale material	metal
Handle material	plastic

The Levenhuk MS3 mechanical scale is designed to move microscope slides on the stage along X axis (left-right) and Y axis (forwards-backwards).

The manufacturer reserves the right to make changes to the product range and specifications without prior notice.

Warranty: six months from date of retail purchase. For further details, please visit our web site: www.levenhuk.com/warranty

The Levenhuk MS3 mechanical scale is designed to move microscope slides on the stage along X axis (left-right) and Y axis (forwards-backwards).

The Levenhuk MS3 mechanical scale is designed to move microscope slides on the stage along X axis (left-right) and Y axis (forwards-backwards).

The Levenhuk MS3 mechanical scale is designed to move microscope slides on the stage along X axis (left-right) and Y axis (forwards-backwards).

The Levenhuk MS3 mechanical scale is designed to move microscope slides on the stage along X axis (left-right) and Y axis (forwards-backwards).

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